

Dear

Thank you again for referring for dental implants. The surgery has been completed successfully and have been advised to return to your office for the restorative procedures.

I am using Implant Direct's ScrewPlant Implants, which come suspended on titanium carriers called "fixture-mounts". These multi-functional, color-coded components are designed for use as implant-level transfers and can also be shortened to create straight abutments. These straight abutments have a groove just above the beveled margin for retention of a snap-on comfort cap and/or plastic transfer.



Click here to view [video](#).

To identify the tooth number associated with each fixture-mount, I have placed the fixture-mount and comfort cap back in the original implant vial and labeled the vial with your name, the patient's name and the tooth number of the implant location. If you use the fixture-mount as a transfer, you can discard the comfort cap.

If you choose to shorten the fixture-mount to create a straight temporary or final abutment, you can use the snap-on comfort cap in the non-esthetic zone or fabricate a temporary crown in the esthetic zone. You would then make a conventional impression or purchase a plastic snap-on transfer for an abutment level transfer. [Link to Straight Abutment Transfer and Analog Flow Chart](#). The video link above demonstrates these procedure options and components.

The labels on the vials, including the one I attached, will provide the laboratory with the information it needs to order the appropriate implant analogs or additional abutments.

I opened an account for you on Implant Direct's Shopping Cart with your address as the "Ship To" and my address as the "Bill To". If you need to order the snap-on transfer and abutment analog for the straight abutment created by shortening the fixture-mount, or if you require a different type of abutment, use this [link to the ScrewPlant abutment ordering page](#). It will take you first to the logon page where you can sign in using your email address and password if you have an account, or you can use the logon page to create an account. At check-out, select "Pre-Order" as the method of payment and insert the patient's name in the order comments so I will be able to approve your charges when the confirmation email is sent to me. Shopping Cart tutorials are available on www.implantdirect.com.

Implant Direct's Customer Service Department can help you with placing an order or to answer any technical questions. Call their toll free number at (888) 649-6425 between 5:00AM and 5:00PM PST Monday through Friday.

If you have any questions regarding this case, please feel free to call me. Thank you again for allowing me to be a part of your Implant Team.

Regards,