



Implant
Direct™

simply smarter.

SX-Plant™ System Catalog



SX-Plant™



4.8mmD
4.1mmD



4.8mmD
4.8mmD



3.7mmD
3.3mmD



6.5mmD
4.8mmD



6.5mmD
5.7mmD



Implant Direct™

simply smarter.

Implant Direct Europe AG
Förlibuckstrasse 150
8005 Zürich, Switzerland
Phone: +41 44 567 8100
www.implantdirect.eu

SX-Plant's All-in-One packaging includes: Cover Screw, Healing Collar, Transfer and Straight Abutment with Snap-on Comfort Cap





"Order Online for High Quality Products at factory-direct prices"



Gerald Niznick, DMD, MSD, Founder and President of Implant Direct

FEATURES AND BENEFITS

1. **Body Design:** Tapered with Mini-threads & progressively deeper, double-lead buttress threads
2. **Surface/Body:** SBM Medium Texture with 17 year history
3. **Surface/Neck:** SBM Light Texture for soft or hard tissue
4. **Diameters:** Standard 4.1mm and 4.8mm plus 3.3mm & 5.7mm
5. **Platforms:** Standard 4.8mm and 6.5mm internal octagon plus 3.7mm hex with Platform Switch
6. **Fixture-mount/Transfer/Abutment:** Engages square within octagon



* The 4.1mm and 4.8mm SX-Plant implants can be inserted into soft bone using Straumann's drills. An additional drill is required for dense bone or countersinking for bone-level placement. Implant Direct's ratchet, insertion tool and 1.25mmD hex tool are also required.