


Implant Direct Pre-August 2007 Sterile Vials with Micro-cracks
Dye Penetration Tests done in 7/07 and Retested 7/09
Study Confirmed that Vials Achieved and Maintained Sterility



TEST 7/07: Hairline crack samples subjected to dye penetration testing in 2007 documented no dye in vial indicating maintenance of sterile seal.

O-Ring maintains sterile seal



Micro-crack does not effect maintainance of sterility



TEST 7/09: Hairline crack samples subjected to dye penetration testing in 2009 documented no dye in vial indicating maintenance of sterile seal.