



**GoDirect One-Piece
Implant System Featuring:**

- One-piece implant in 3 diameters with Zest Locator® compatible platform
- Tapered, self-tapping body for increased stability in soft bone
- 3.0mm "mini-implant" with 1.8mm apical end allowing insertion with 2.3mmD Drill
- Double-lead threads for faster insertion, with mini-threads for reduced stress concentration near crest of bone
- Internal threads, capped off with small internal cover-screw, allows for conversion to multi-unit abutment for fixed-detachable bars or bridges
- All-in-one packaging includes snap-on transfer and comfort cap





GoDirect


GoDirect: One-Piece Implant for One-Stage Surgery

SBM Surface Implants

Includes Extender, Carrier/Transfer and Comfort Cap


3.0mmD

Code	Diameter	Length	Platform	Price
423010	3.0mmD	10mmL	4.0mmD	\$150.00
423011	3.0mmD	11.5mmL	4.0mmD	\$150.00
423013	3.0mmD	13mmL	4.0mmD	\$150.00
433010	3.0mmD	10mmL	4.0mmD	\$150.00
433011	3.0mmD	11.5mmL	4.0mmD	\$150.00
433013	3.0mmD	13mmL	4.0mmD	\$150.00



3.7mmD

Code	Diameter	Length	Platform	Price
423710	3.7mmD	10mmL	4.0mmD	\$150.00
423711	3.7mmD	11.5mmL	4.0mmD	\$150.00
423713	3.7mmD	13mmL	4.0mmD	\$150.00
433710	3.7mmD	10mmL	4.0mmD	\$150.00
433711	3.7mmD	11.5mmL	4.0mmD	\$150.00
433713	3.7mmD	13mmL	4.0mmD	\$150.00



4.7mmD

Code	Diameter	Length	Platform	Price
424708	4.7mmD	8mmL	4.0mmD	\$150.00
424710	4.7mmD	10mmL	4.0mmD	\$150.00
424713	4.7mmD	13mmL	4.0mmD	\$150.00



Locator Abutment Components

Code	Retention Force	Description	Price
1000-92C	5.0lbs	Replacement Male Clear - 4ea	\$21.00
1000-92P	3.0lbs	Light Replacement Male Pink - 4ea	\$21.00
1000-92B	1.5lbs	Extra Light Replacement Male Blue - 4ea	\$21.00
1000-92G	4.0lbs	Extended Range Male Green - 4ea	\$21.00
1000-92R	1.0lb	Extra Light Extended Range Male Red - 4ea	\$21.00
1000-93	-	Black Processing Male - 4ea	\$21.00
1000-98	Assorted	Abutment Processing Cap/Male Kit - 2ea	\$58.00



Locator Instrumentation

Code	Description	Price
LDT	Locator Abutment Driver Tool - Long	\$44.00
LDTs	Locator Abutment Driver Tool - Short	\$44.00
LDD	Locator Abutment Driver Drill - Long	\$44.00
LDDs	Locator Abutment Driver Drill - Short	\$44.00
1000-96	Core Tool - Male Removal & Seating Tool and Abutment Driver	\$75.00

