



ScrewIndirect One-Piece Implant System Featuring:

- **Screw Receiving head with retentive groove & beveled shoulder 2.5mm above bone**
- **Double-lead threads for 2x faster insertion and quadruple-lead threads near top**
- **Color coded plastic fixture mount usable as a snap-on transfer**
- **Snap-on plastic comfort cap**
- **Common prosthetic platform across all four diameters of implants**
- **Includes 2mm implant extender**
- **Tapered body allows insertion using Zimmer Dental's Tapered Screw-Vent® drills**

5.0mmD Platform



5.0mmD Platform



5.0mmD Platform




5.0mmD Platform




3.0mmD Implant




3.7mmD Implant



4.7mmD Implant





5.7mmD Implant




ScrewIndirect Packaging

6.0mmD Platform
403713
13mm LONG

Carrier/Transfer



Extender Transfer



Snap-On Comfort Cap



Extender

Snap-On Comfort Cap and Transfer

Seat retentive snap-on comfort caps




Seat retentive snap-on transfers




Clinical Applications

Screw-Retained Bar Overdenture



Fixed-Detachable Restoration





ScrewIndirect: One-Piece Implant for One-Stage Surgery

Indirect

SBM Surface Implants

Includes Extender, Carrier/Transfer and Comfort Cap

3.0mmD

Code	Diameter	Length	Platform	Price
403008	3.0mmD	8mmL	5.0mmD	\$150.00
403010	3.0mmD	10mmL	5.0mmD	\$150.00
403013	3.0mmD	13mmL	5.0mmD	\$150.00
403016	3.0mmD	16mmL	5.0mmD	\$150.00




3.7mmD

Code	Diameter	Length	Platform	Price
403708	3.7mmD	8mmL	5.0mmD	\$150.00
403710	3.7mmD	10mmL	5.0mmD	\$150.00
403713	3.7mmD	13mmL	5.0mmD	\$150.00
403716	3.7mmD	16mmL	5.0mmD	\$150.00




4.7mmD

Code	Diameter	Length	Platform	Price
404708	4.7mmD	8mmL	5.0mmD	\$150.00
404710	4.7mmD	10mmL	5.0mmD	\$150.00
404713	4.7mmD	13mmL	5.0mmD	\$150.00
404716	4.7mmD	16mmL	5.0mmD	\$150.00



5.7mmD

Code	Diameter	Length	Platform	Price
405708	5.7mmD	8mmL	5.0mmD	\$150.00
405710	5.7mmD	10mmL	5.0mmD	\$150.00
405713	5.7mmD	13mmL	5.0mmD	\$150.00
405716	5.7mmD	16mmL	5.0mmD	\$150.00



Components for use with 2mmL ScrewIndirect Extender on pg 22

ScrewIndirect Implant Copings

Includes Fixation Screw [1000-01]

Code	Platform	Description	Price
4000-64	5.0mmD	Gold Coping	\$75.00
4000-65	5.0mmD	Titanium Coping	\$40.00
4000-66	5.0mmD	Plastic Coping	\$35.00
4000-67	5.0mmD	Gold & Plastic Coping	\$75.00
4000-69	5.0mmD	Titanium Temporary Coping	\$40.00
1000-01	-	Replacement Fixation Screw - 2ea	\$30.00



ScrewIndirect Comfort Cap, Transfers & Analogs

Code	Platform	Description	Price
4000-13	5.0mmD	Plastic Comfort Cap - 2ea	\$30.00
4000-02	5.0mmD	Plastic Implant Transfer - 2ea	\$30.00
4000-62	5.0mmD	Open-Tray Transfer plus Fixation Screw	\$30.00
4000-63	5.0mmD	Closed-Tray Transfer plus Fixation Screw	\$30.00
4000-06	5.0mmD	Aluminum Implant Analog	\$20.00
1000-10	-	Replacement Open-Tray Fixation Screw - 2ea	\$30.00
1000-16	-	Replacement Closed-Tray Fixation Screw - 2ea	\$30.00

