



**Implant
Direct™**

simply smarter.

Implant Direct Int'l
8840 W. Russell Road #210,
Las Vegas, NV 89148
Manufacturing/Sales: 1-818-444-3333
www.implantdirect.com



**“Order Online for High Quality Products
at Factory-Direct Prices”**

Gerald Niznick, DMD,MSD, Founder and President of Implant Direct

GoDirect™ (Pat. Pend.)


**One-Piece Implant with
*Locator® Compatible Platform**
All-in-One packaging includes
Snap-on Transfer and Comfort Cap
USA List Price = \$150





1-PIECE IMPLANT DESIGN FEATURES

1. Platform: Matches Zest Locator® Abutment
2. Body: Tapered with self-tapping cutting grooves
3. Surface: Medium Rough Texture - SBM 17 year history
4. Lengths: 10mm, 11.5mm & 13mm
5. Diameters: 3.0mmD, 3.7mmD & 4.7mmD
6. Collar Height: 1.5mm & 3.0mm


Snap-on
Comfort Cap
Included



Snap-on
Transfer
Included



Zest Locator®
Male Attachment
Not Included



*Locator® is a registered Trademark of Zest Anchor Company