

Table of Contents



Legacy™ System
page 2-5



Compatible with
*Zimmer Screw-Vent®




Legacy3 3.2/3.0 Legacy3 3.7/3.5 Legacy3 4.2/3.5 Legacy3 4.7/4.5 Legacy3 5.2/4.5 Legacy3 5.7/5.7

Legacy3 Implants with progressively deeper threads.
Packaging includes cover screw, healing collar, transfer and abutment



Spectra-System®
page 6-9




ScrewPlant® ScrewPlus® ScrewDirect® ScrewIndirect® GoDirect™ ScrewRedirect®

Application specific 1- and 2-piece Implants with all-in-one packaging
Same implant body and surgical protocol for all six Implants



Tri-Lobe System
page 10-13



Compatible with
*Nobel Biocare™



ReActive® 3.7/3.5 ReActive 4.2/3.5 ReActive 4.7/4.3 ReActive 5.7/5.0

ReActive Implants with progressively deeper threads.
Packaging includes cover screw, transfer and abutment



SwissPlant® System
page 14-17



†Compatible with
*Straumann®



SwissPlant 3.3/3.7 SwissPlant 4.1/4.8 SwissPlant 4.8/4.8 SwissPlant 4.8/6.5 SwissPlant 5.7/6.5

SwissPlant Implants with progressively deeper threads.
Packaging includes cover screw, healing collar, transfer and abutment