











































Tool Application

Tool Used	Compatible Parts Used With Specified Tool											
Ratchet (R)			Hex Tool 2.0 (L & S)		Hex Tool 2.5 (L & S)		Hex Tool 3.0 (L & S)		All Fixture Mounts: ScrewPlant/Plus, RePlus & Legacy			
												
Stainless Steel Handle (SSH)			Hex Tool 2.0 (L & S)		Hex Tool 2.5 (L & S)		Hex Tool 3.0 (L & S)		All Fixture Mounts: ScrewPlant/Plus, RePlus & Legacy			
												
Fixture Mount Drill (FMD)			Hex Tool 2.0 (L & S)		Hex Tool 2.5 (L & S)		Hex Tool 3.0 (L & S)		All Fixture Mounts: ScrewPlant/Plus, RePlus & Legacy			
												
Torque Wrench (TW-30)	Hex Tool 1.25											
												
Torque Wrench (TW-35)									All Fixture Mounts: ScrewPlant/Plus, RePlus & Legacy			
												



27030 Malibu Hills Road
 Calabasas Hills, CA 91301

Phone: 818.444.3333 opt.2 Fax: 818.444.3400

t e c h @ i m p l a n t d i r e c t . c o m