




























Hex Tool Application

Hex Tool Used		Compatible Parts Used With Specified Hex Tool				
1.25 Hex Tool		All Cover Screws	All Healing Collars	All Ball Abutments	All Screw-Receiving Abutments	All Abutment Screws
	 Hex Tool (L or S)					
2.0 Hex Tool			ScrewDirect 3030 & 3037 Implants			All ScrewRedirect Implants
	 Hex Tool (L or S)		 			  
2.5 Hex Tool		All ScrewPlant/Plus FMT's, 1037/ 47, & 2030/ 37 Implant Diameters	ScrewDirect 3047 & 3057 Implants	Legacy 8037 & 8047 Implants	All RePlus FMT's	All ScrewIndirect Implants
	 Hex Tool (L or S)	   	 	 	  	  
3.0 Hex Tool		ScrewPlant/Plus Implants 1057, 2047 & 2057		Legacy 8057 Implants		
	 Hex Tool (L or S)	  				



27030 Malibu Hills Road
Calabasas Hills, CA 91301

Phone: 818.444.3333 opt.2 Fax: 818.444.3400

tech@implantdirect.com