


Subject:
Product Release: Update to Fixture Mount Drill

Date: 9-12-08




To: Implant Direct Customers



The updated version of the Fixture Mount Drill is now available!!

Implant Direct introduces the new and improved version of our Fixture Mount Drill [FMD]. The drill is designed with a sterilizable surgical plastic split O-ring at its apical end. This plastic split O-ring replaces the current rubber O-ring. During insertion of the Fixture Mount Drill over the mating square connection of the fixture mount, the O-ring will expand and create a frictional contact with the outer walls of the fixture mount carrying the implant. This frictional connection between the two mating surfaces allows the clinician to carry the corresponding implant assembly to the surgical site while connected to the surgical handpiece. Can be used in combination with the Drill Extender [DE].

The USA list price for the Fixture Mount Drill is \$80.00.



27030 Malibu Hills Road, Calabasas Hills, CA 91301

Phone: 818.444.3333 opt2: 888-NIZNICK (649-6425) Fax: 818.444.3400

w w w . i m p l a n t d i r e c t . c o m