

Sun, 02/21/2010 - 16:34 — [Tom House DDS MD](#)

[Nobel Implant System.... No reason to use another](#)

I have been using Nobel for years. Started with the Mark IV's now I use the Nobel Replace tapered and the Nobel active. Surface is great. Very nice healing. My personal success rate is over 95%. Restoratively it is easy to restore.

Mon, 02/22/2010 - 01:52 — [niznick](#) [reply](#) [Here are a few reasons:](#)

1. \$350 per implant plus \$48 for cover Screw. You can get an Implant Direct RePlant that goes in with the same Nobel drills, accepts the same abutments, has the same color coded platform, and comes with a cover screw and fixture-mount that is a transfer and temporary abutment - \$150

[http://www.implantdirect.com/newsletter\\_0/newsletter\\_18/ReActive\\_ad.pdf](http://www.implantdirect.com/newsletter_0/newsletter_18/ReActive_ad.pdf)

2. RePlant is made of titanium alloy so the 3.5mmD implant will not break during insertion as has been reported with the Nobel Replace, made of CP titanium.

3. RePlant is self-tapping so you do not need to use a bone tap, saving a step and assuring increases stability.

[http://www.implantdirect.com/us/pop-out/RePlant\\_Surgical&Procedures\\_HiQ.htm](http://www.implantdirect.com/us/pop-out/RePlant_Surgical&Procedures_HiQ.htm)

4. RePlant has micro-threads that add to the surface area and to initial stability, while reducing bone loss.

[http://www.implantdirect.com/newsletter\\_0/newsletter\\_18/newsletter\\_20\\_tg\\_a.htm#patent](http://www.implantdirect.com/newsletter_0/newsletter_18/newsletter_20_tg_a.htm#patent)

5. RePlant showed only 1% failure compared to 13% with Nobel Replace in a controlled clinical study where both implants were placed side-by-side in the same edentulous jaws followed by immediate loading. - [http://www.implantdirect.com/pdf/EAO\\_ImplGuidSurg\\_Final.pdf](http://www.implantdirect.com/pdf/EAO_ImplGuidSurg_Final.pdf)

6. RePlant implants from Implant Direct have demonstrated that they are packaged in vials that maintain seal integrity, assuring that the implants remain sterile for the entire 5 years claimed on the label – Nobel Biocare's Implants have repeatedly failed such standard dye leak tests.

[http://www.implantdirect.com/newsletter\\_0/newsletter\\_18/nobel.htm](http://www.implantdirect.com/newsletter_0/newsletter_18/nobel.htm)

7. An FDA Medwatch Report filed by a dentist documented 38.6% failures compared to only 3.8% failure with Implant Direct's RePlant Implants.

[http://www.implantdirect.com/newsletter\\_0/newsletter\\_18/nobel.htm#medwatch2](http://www.implantdirect.com/newsletter_0/newsletter_18/nobel.htm#medwatch2)

8. Implant Direct's RePlus and ReActive Implants offer All-in-1 Packaging that includes the cover screw, transfer and final abutment for about 70% less than similar parts from Nobel Biocare.

[http://www.implantdirect.com/us/pop-out/RePlus\\_Packaging\\_HiQ.htm](http://www.implantdirect.com/us/pop-out/RePlus_Packaging_HiQ.htm)

[http://www.implantdirect.com/us/pop-out/ReActive\\_Packaging\\_HiQ.htm](http://www.implantdirect.com/us/pop-out/ReActive_Packaging_HiQ.htm)

9. Implant Direct's Nobel Replace compatible abutments and drills cost 50-70% less

10. Implant Direct's RePlus and ReActive tri-lobe implants with final abutments included have platform switching between the top of the implant and the tri-lobe abutment.

[http://www.implantdirect.com/us/pop-out/ReActive\\_RestorativeProcedures\\_FixtureMount\\_HiQ.htm](http://www.implantdirect.com/us/pop-out/ReActive_RestorativeProcedures_FixtureMount_HiQ.htm)