





**Implant Direct**  
simply smarter.  
Sybron Dental Specialties

# RePlant® System Catalog



ReActive®



RePlant®  
Internal Tri-Lobe



RePlus®  
Internal Tri-Lobe



ReActive™  
Internal Tri-Lobe



www.implantdirect.com  
 Implant Direct Sybron International  
 1-818-444-3333 • 1-888-NIZNICK



**"Order Online for High Quality Products at Factory-Direct Prices"**

Gerald Niznick, DMD, MSD, Founder and President

### ReActive™ Tri-Lobe Implant

Platform Shifting, Micro-threads", Deeper Threads and Round Apex.

All-in-one Packaging Includes: Cover Screw, Transfer and Straight Abutment.

ReActive USA price = \$200

NobelActive USA price = \$694



Three Implant Designs and Packaging Options  
 Prosthetically Compatible to NobelReplace Implant  
 Implants include Cover Screw

Select Implant based on Price, Packaging and Thread Design

1. RePlant: \$150 includes Temporary Abutment/Transfer
2. RePlus: \$175 includes "Snappy" Abutment/Transfer
3. ReActive: \$200 includes Preparable Abutment/Transfer

**RePlant® \***  
 Body Dimensions  
 Matches Replace



**RePlus®**  
 Body Taper  
 For Bone Expansion



**ReActive™**  
 Buttress Threads  
 For Increased Surface




\*RePlant implant surgically compatible with NobelReplace® Drills

#### Free Implants for Nobel Biocare Customers


3 Surgically Compatible, RePlant® Implants with Titanium Transfer/Temporary Abutments plus 2 Contoured Abutments

US List Price for 3 Implants, Abutments and Prosthetic Components:

**RePlant**  
 = \$620  
**Nobel Biocare®**  
 = \$1922



FREE Introductory Offers at [www.implantdirect.com](http://www.implantdirect.com)  
 Terms and Conditions Apply



Compare US List Price	Laboratory Abutment	Straight Snap-On	Straight Contoured	15° Angled Contoured	Gold/Plastic	Zirconia Abutment	Plastic Temporary Abutment	Ball Attachment	GPS™ Attachment	GPS™ Angled Attachment	Multiple-Unit w / Cap & Transfer	Multiple-Unit w / Cap & Transfer Angled
Implant Direct	\$75	\$85	\$85	\$85	\$100	\$120	\$35	\$107	\$100	\$120	\$85	\$100
Nobel Biocare	\$145	\$186	\$170	\$225	\$216	\$285	\$72	\$217	\$149	N/A	\$229	\$274
Cost Savings	48%	54%	50%	62%	54%	58%	51%	51%	33%	N/A	63%	64%

