




Legacy Prosthetic Vial and Cap Label



- | | |
|--|---|
| | 1) Implant direct products carry the CE mark and fulfill the requirements of the Medical Device Directive 93/42/EEC |
| | 2) Product code or reference number |
| | 3) Product description |
| | 4) Use before expiration date (XXXXXXX) |
| | 5) Gamma sterilization symbol |
| | 6) Single-use only symbol |
| | 7) Caution: US law restricts the sale of these devices by or on the order of a dentist or physician |
| | 8) Refer to Instructions For Use (IFU) |
| | 9) Lot number (XXXXXXXXXX) |
| | 10) Implant body diameter color reference on surgical tray |
| | 11) Implant prosthetic platform diameter (written and color coded) |
| | 12) Non-sterile |