




Overdenture Abutments Now Available for Legacy3 Narrow-Diameter Implants

The versatility of the 3.2mm-diameter Legacy3 implant, which features a 3.0mm platform, has just increased thanks to the availability of overdenture abutments.

Narrow-diameter implants are a viable alternative to traditional implants in edentulous cases. Narrow-diameter implants may be indicated when the mandibular ridge is limited, the patient prefers conservative treatment or other health factors prevent extensive surgical procedures.

Effective immediately, overdenture abutments for the 3.2mmD Legacy3 implant are available in three different collar heights. This range allows you to select the best fit depending upon the amount of soft tissue present in each case. Please note that these overdenture abutments are recommended for use only in multi-unit cases.

Make the most of narrow-diameter implants with Legacy3's many abutment choices!



3.0mmD platform Legacy3 Overdenture Abutments

8530-61, 8530-62, 8530-64

Expand your treatment options now!

[Watch Legacy3 Online Tutorial](#)

[Shop Online Now](#)