

Subject:
Product Release: Dual-Grip™ Narrow Healing Collars

Date: 1-26-09


To: Implant Direct Customers

Narrow Diameter Color Coded Healing Collars for Replant System

Implant Direct introduces the **Narrow Diameter Healing Collars** manufactured with the unique **Dual-Grip™** configuration. This feature enables clinicians to use either the standard 1.25mmD Hex Tool OR the Nobel BioCare™ compatible **UniGrip™** driver.

The Healing Collars are manufactured from Titanium Alloy and are color coded to match the colored platforms of the Replant, RePlus, ReActive and Nobel BioCare™ tri-lobe interface systems.

These narrow diameter components follow the same outside diameter as the interface to which they connect, namely; 3.5mmD, 4.3mmD, 5.0mmD and 6.0mmD. Available in two lengths, 3mmL and 5mmL corresponding to the currently available lengths of other healing components in the system.



27030 Malibu Hills Road, Calabasas Hills, CA 91301

Phone: 818.444.3333 opt2: 888-NIZNICK (649-6425) Fax: 818.444.3400


w w w . i m p l a n t d i r e c t . c o m

Standard Healing Collars- Require 1.25mmD Hex Driver

Healing Collars

Compatible to use with RePlant, RePlus and ReActive Implants

Code	Width	Length	Platform	Price
6035-13	4.6mmD	3mmL	3.5mmD	\$30.00
6035-15	4.6mmD	5mmL	3.5mmD	\$30.00
6043-13	5.6mmD	3mmL	4.3mmD	\$30.00
6043-15	5.6mmD	5mmL	4.3mmD	\$30.00
6050-13	6.1mmD	3mmL	5.0mmD	\$30.00
6050-15	6.1mmD	5mmL	5.0mmD	\$30.00
6060-13	6.6mmD	3mmL	6.0mmD	\$30.00
6060-15	6.6mmD	5mmL	6.0mmD	\$30.00




Dual-Grip™ Narrow Healing Collars- Use either 1.25mmD Hex or UniGrip Driver

Narrow Healing Collars

Components include Dual-Grip™ feature for use with 1.25mmD Hex Tool or Nobel Unigrip ScrewDriver



Code	Width	Length	Platform	Price
6035-13N	3.5mmD	3mmL	3.5mmD	\$30.00
6035-15N	3.5mmD	5mmL	3.5mmD	\$30.00
6043-13N	4.3mmD	3mmL	4.3mmD	\$30.00
6043-15N	4.3mmD	5mmL	4.3mmD	\$30.00
6050-13N	5.0mmD	3mmL	5.0mmD	\$30.00
6050-15N	5.0mmD	5mmL	5.0mmD	\$30.00
6060-13N	6.0mmD	3mmL	6.0mmD	\$30.00
6060-15N	6.0mmD	5mmL	6.0mmD	\$30.00



Dual-Grip™ Healing Collars require 1.25mmD Hex Driver **or** Industry Standard UniGrip™ Driver