

**Subject:**  
**Product Release: RePlus® HA**

**Date: 4-09-08**

**To: Implant Direct Customers**

**RePlus HA coated implants are now available!!!!**

Produced in three standard implant diameters: 3.7mm, 4.7mm and 5.7mm, the implant body tapers evenly from top to bottom to facilitate insertion into an undersized osteotomy, in soft bone this provides compression that assists in increasing initial stability. RePlus implants are placed using the standard Spectra System drill kit from Implant Direct.

Manufactured from high strength titanium alloy, the implants tri-lobe color-coded platforms, magenta (3.5mmD), gold (4.3mmD) and blue (5.0mmD) are compatible with Nobel Replace® abutments. Implants have a microtextured surface from the top of the implant excluding the color coded prosthetic interface. This surface treatment transitions into a high crystalline HA coating below the mini threads (*3mm from the top of implant*) extending to the apex of the implant. The high quality, cost effective RePlus HA implants have all the design features and advantages of the well received RePlus SBM, combined with the bioactive advantages of a HA surface coating.

**The USA list price for the RePlus HA is \$225.00.**



3.7mmD Body  
713708 - 713716  
3.5mmD Platform



4.7mmD Body  
714708 - 714716  
4.3mmD Platform



5.7mmD Body  
715708 - 715716  
5.0mmD Platform



27030 Malibu Hills Road, Calabasas Hills, CA 91301

Phone: 818.444.333 opt2: 888-NIZNICK (649-6425) Fax: 818.444.3400

**w w w . i m p l a n t d i r e c t . c o m**